



# PANASIL<sup>®</sup> MONOPHASE MEDIUM



**Panasil monophase Medium** is a medium-viscosity monophase impression material based on A silicone. It is available in the comfortable 5:1 jumbo cartridge as well as the familiar 50 mL cartridge and is characterized by a particularly high initial hydrophilicity for precision in extreme situations.



## Precise impressions

- Flows optimally and, thanks to its exceptional initial hydrophilicity, also ensures reliability and precision in difficult oral situations.
- Dimensionally stable elastic recovery properties.

## Universally applicable

- Just as suitable for preparing crowns/bridges, inlays/onlays, and veneers as for fixation impressions.

## Perfectly customized


- The practical, short intraoral setting time (2 minutes) combined with high final hardness (Shore A 60) enable particularly precise impressions of the three-dimensional relationships in the fixation.





## Panasil monophase Medium medium viscosity

- + Monophase impressions
- + Fixation impressions
- + Functional impressions
- + Pick-up impressions
- + Double-mix impressions
- + Reline impressions


- + Highly recommended
- + Recommended




- |                                      |   |       |
|--------------------------------------|---|-------|
| Dynamic mixers                       |  | p. 54 |
| Mixing tips, green, 100 tips         |   | p. 54 |
| Electrical dosing and mixing devices |   | p. 56 |
| Panasil <sup>®</sup> Adhesive        |   | p. 53 |

380 mL **Intro pack 5:1**   
380 mL cartridge,  
10 dynamic mixers

REF 14708

760 mL **Refill pack 5:1**   
2 x 380 mL cartridges

REF 14709

100 mL **Normal pack 1:1**   
2 x 50 mL cartridges, 6 green mixing tips

REF 13501